
97050164
rev. 011
06/2014

ITALIANO	3
ENGLISH.....	25
FRANÇAIS	46
DEUTSCH	69
ESPAÑOL.....	93
PORTUGUÊS.....	117
РУССКИЙ	139
POLSKI.....	163
中文.....	185

ZEN-X / ZEN-XI

IT
EN
FR
DE
ES
PT
RU
PL
ZH

INDICE

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CONTENTS

1 Foreword.....	26
1.1 Conventions in text	26
2 General information.....	26
2.1 Warnings for use.....	27
2.2 Safety and Health Care	28
2.3 Maintenance and Disposal	28
3 Standards	29
4 Package Contents.....	29
5 Product Identification.....	30
6 Warranty	33
6.1 Software not covered by warranty	33
6.2 Limitation of liability.....	33
7 Minimum system requirements.....	34
8 Technical characteristics	34
9 Sensor characteristics	35
9.1 Compatibility with x-ray generators.....	35
10 Disposable covers	37
11 Cleaning and disinfecting	37
12 Use of the centering device	39
13 Installing the sensor.....	39
14 Use of the sensor.....	39
15 Turning the sensor on and off	40
16 Acquiring an x-ray image	40
17 Status light	41
18 ZEN-X.....	42
19 ZEN-Xi	42
20 Quality of the x-ray images	44
21 Problem solving	44

1 Foreword

The ZEN-X / ZEN-Xi digital system has been developed to simplify the entire intraoral x-ray acquisition procedure and display the images on a computer screen. Thanks to the new ergonomic design, MyRay® sensors permit easy intraoral positioning. The smoothed edges and rounded corners comfortably adapt to the shape of the patient's mouth, allowing them to be easily positioned. MyRay sensors are available in two interchangeable sizes to satisfy various diagnostic requirements.

The ZEN-X / ZEN-Xi electronic module is compatible with the extremely fast USB 2.0 standard, thereby cutting down the time that elapses between x-ray exposure and display of the image on the computer screen to a few seconds.

The product is available in two versions that adopt the same technology but differ for the appropriate application:

- Zen-X is designed as a portable device;
- Zen-Xi is constructed to be integrated in dental units that are mechanically compatible.

In Zen-X, the USB connection makes the system convenient and portable. In fact, no power supply adapters are needed as power is fed directly through the USB port thanks to the low consumption requirements. The handheld control can be easily moved from one room to another and from a desktop computer to a notebook.

In Zen-Xi, the x-ray sensor is placed on-board the dental unit like any other instrument on the instrument board. S

A computer and program to view x-ray images are required for use of both versions. If used in conjunction with a dental office practice management software, the x-ray images can be associated to each patient and saved for processing and viewing when required.

The ZEN-X / ZEN-Xi system utilizes a communication standard called TWAIN®, adopted in numerous electronic devices such as scanners and digital cameras. TWAIN® guarantees compatibility with all the best programs for managing, processing and storing digital images.

Regardless of the selected program, whether it is iRYS produced by MyRay or another third party program, refer to the manual supplied together with the program for all the warnings, precautions and operating instructions.

ZEN-X / ZEN-Xi is supplied together with software called iCapture, that assures the x-ray images are correctly transferred from the electronic module to the computer.

1.1 Conventions in text

This User's manual refers to both the Zen-X and Zen-Xi systems, providing specific details only when there are differences between the two versions.

The manual refers to a "computer", Personal Computer, "work post", WorkStation or WS indifferently.

2 General information

Pay special attention to the sections of the manual in which the indicated symbols appear:

	Warnings regarding operator or patient safety
	Warnings regarding risk of product damage or failure to operate, important warnings regarding the warranty.

The ZEN-X / ZEN-Xi system and relative iCapture software are developed and produced by Cefla S.C. - Cefla Dental Group, Imola (IT), herein referred to as the Manufacturer.

	The instructions inform the user on how to properly operate the ZEN-X / ZEN-Xi system. As regards the iCapture software instructions, consult the specific manual. Carefully read both manuals before attempting to use the sensor and program.
--	---

In order to use the ZEN-X / ZEN-Xi system, software for capturing and saving the images is needed that is not part of the ZEN-X / ZEN-Xi system. Consult the relative manual for information about installation and use of the image management software.

It is forbidden to reproduce, store and transmit this publication by any means (electronically, mechanically, by photocopying, translating or in other ways) without a written authorisation issued by the Manufacturer.

Any information, technical specifications and illustrations contained in this manual are not binding on the Manufacturer. The Manufacturer reserves the right to introduce modifications and technical improvements without having to modify these instructions.

Important: in accordance with privacy laws in force in several countries, all sensitive personal information must be adequately protected. In addition, patients must sign a consent form before personal information or images are transmitted across networks.

The original text of this manual is in Italian.
Translation from the original in Italian

All the registered trade marks and product names mentioned in this manual are the property of their respective owners.

2.1 Warnings for use

The ZEN-X / ZEN-Xi system is designed to work only when connected to accessories equipped with adequate interface software. For this reason, neither the sensor, electronic interfaces nor software components ("drivers" installed in computer and "firmware" in the devices) are compatible with other commercial devices. Therefore use of the ZEN-X / ZEN-Xi system and relative software in conjunction with other commercial devices is not assured or recommended

The ZEN-X / ZEN-Xi system uses the TWAIN® protocol for data transmission. It can be used with any program able to acquire images from TWAIN® devices (e.g. scanner – digital cameras). Medical programs should be used as they assure data security and quality of the images.

Even though other software interfaces are compatible with the ZEN-X / ZEN-Xi system and relative software components, it is not advisable to use other x-ray image acquisition software simultaneously in the same computer employed to capture images with the ZEN-X / ZEN-Xi system or use other software simultaneously for acquiring images in general (scanners, digital cameras, etc).

Some manufacturers of dental office management programs protect their products by purposely making them incompatible with equipment manufactured by third parties. For this reason, it is not possible to guarantee full compatibility of the ZEN-X / ZEN-Xi system with all programs currently available.

We recommend to regularly make backup copies of all acquired images.

The PC should have adequate antivirus software and be used only as a work instrument.

Installing new programs in the computer or updating the operating system may interfere with the TWAIN® driver or image capturing software. After installing new programs in the computer or updating the operating system check system operation before attempting to use on a patient.

Electronic apparatuses may cause or be subject to interferences when used near other electromagnetic equipment such as mobile phones, personal computers equipped with Wireless LAN cards and microwave ovens. Keep the parts of the ZEN-X / ZEN-Xi system and Personal Computer used to acquire and save the x-ray images away from RF sources such as wireless LAN cards, other radiation devices, Home RF devices, microwave ovens; the recommended distance is at least 1 meter, 2 meters in the case of microwave ovens.

Important: in the event of Personal Computer failure while the x-ray image is being transferred (software "crash"), in many cases the x-ray picture is stored in the electronic interface until it is successfully transferred or the interface is shut off or disconnected. The manual procedure described in paragraph "Recovering the last image acquired" (iCapture user instructions) can be carried out to recover the image. This event is extremely unlikely as it takes just a few seconds to transfer the image from the electronic control to the computer.

Consult the manuals that deal with the individual apparatuses for instructions on how to use the equipment in conjunction with the ZEN-X / ZEN-Xi system (computer, x-ray unit, etc).

Only specially trained technicians should install the other system components (computer or computer network, soft-

ware to manage and save images, x-ray generator, etc...). In particular, keep in mind that installation of the x-ray equipment is checked and inspected by a qualified technician.

	The USB® connection required for ZEN-X / ZEN-Xi system operation is not a normal electric connection. In fact, special cables are needed (marked USB® HiSpeed®, see paragraph “Product identification”). To assure trouble-free operation, a single USB® 2.0 cable cannot be longer than 4.5m. If the ZEN-X / ZEN-Xi system needs to be installed with longer cables a USB® repeater cable must be used for each 4.5m section (maximum) and for at most three sections (two repeaters).
--	---

	The sensor used to capture images is fragile and bothered by static electricity. Handle it with care. Do not deform or squeeze it. Do not touch the electric contacts when the connector is not plugged into the electronic control module. Do not disconnect it while the interface is on; see paragraph “Using the sensor”.
--	--

2.2 Safety and Health Care

- The ZEN-X / ZEN-Xi digital system is a medical device for intraoral x-ray acquisition. It is intended for use only by qualified dental professionals. Do not use the system for any purposes other than capturing intraoral radiographs. The device can be used only by clinicians well instructed in the fields of dentistry and radiology.
- Never use electronic apparatuses in the vicinity of life support (e.g: pacemaker or cardiac stimulators) and hearing aids. In medical facilities, before using any electronic devices, always assure they are compatible with the equipment present.
- **The disposable protection control sheaths provided are mandatory to safeguard the patients against cross-contamination. The disposable protection control sheaths are class IIa medical devices. They cannot be replaced with others that are poorer quality. Contact the dealer who supplied the sensor to obtain additional disposable protection control sheaths.**
- Cover all the components which personnel may touch with disposable protection as they may be contaminated by indirect contact with the patient's mouth. In particular, be careful when handling the Personal Computer mouse, keyboard or touch screen.
- Never use in the presence of mixtures of flammable anesthetic gas with air, oxygen or nitrogen protoxide.
- Some parts (USB® cable, silicon protection, disposable infection control sheaths, parts of the positioner, items used for packaging, x-ray sensor) may result in suffocation if swallowed or used erroneously and inappropriately. Avoid unintended, inappropriate and misuse and keep out of reach of children.

2.3 Maintenance and Disposal

The equipment does not contain any parts that can be repaired directly by the user. In the event of malfunctions, do not try to perform any maintenance jobs. Contact the manufacturer or local dealer directly at the phone numbers given on the warranty certificate. If the apparatus has to be returned to the manufacturer or Service Center for any reason, completely disinfect the outside of the apparatus with a specific product (see paragraph “Cleaning and disinfecting”) and send it back preferably in its original box.

No electronic parts of the ZEN-X / ZEN-Xi system require maintenance. Devices may be broken, the protective means for electric safety are violated and the warranty is void if the sensor or interface casings are opened to reach the circuits inside.

Do not use the sensor on a patient if a system malfunction is present or suspected.

At the end of the appliance working life, dispose of in accordance with the applicable regulations in force. We recommend to disinfect all the appliance external parts before disposal and to sort out the stripped materials for separate waste collection.

In compliance with Directives 2002/95/EC, 2002/96/EC and 2003/108/EC regarding restriction of the use of certain hazardous substances in electrical and electronic equipment along with waste electrical and electronic equipment, it is forbidden to dispose of this equipment in the municipal waste stream as unsorted municipal waste. This equipment must be sorted and disposed of separately. When new equipment that is similar is purchased, the old equipment must be given to the dealer for disposal. As regards reuse, recycling and other forms of recovery of waste electrical and electronic equipment, the producer carries out the functions set forth by laws and regulations currently in force in the individual countries. A high level of separate collection of waste electrical and electronic equipment is indispensable to efficiently recycle, treat and dispose of the equipment. Recycling and treatment operations should comply with minimum standards to assure human health and high environmental protection as well as favor recycling of the materials included in the equipment. The symbol indicating separate collection for electrical and electronic equipment consists of the crossed out bin marked on the equipment. Illegal waste clearance and disposal shall be punished as established by laws and regulations currently in force in the individual countries.

3 Standards

The ZEN-X / ZEN-Xi system is a class IIa medical device for taking intraoral x-rays.

The ZEN-X / ZEN-Xi system and its accessories are designed and constructed in compliance with the following standards:

93/42/CEE and subsequent changes	Electromedical appliance directive
EN 60601-1	General regulations for electromedical appliance safety
EN 60601-1-2	Electromedical appliances: Electromagnetic compatibility – Prescriptions and tests
EN 60878	Graphic symbols for electromedical appliance functions
UNI CEI EN ISO 14971	Medical appliances: Risk management application to medical appliances


4 Package Contents

The package contains:

ZEN-X package (stand alone)


	ZEN-X system
	CD-Rom with the software and drivers
	Warranty certificate
	Compliance certificate

ZEN-Xi package (in dental unit)

	Interface for instrument board with sensor cradle and cable ring
	CD-Rom with the software and drivers
	Warranty certificate
	Compliance certificate

Sensor package

	X-ray sensor
	MyRay® positioners introductory kit including package of disposable infection control sheaths
	Sensor warranty certificate

5 Product Identification
Intraoral sensor

The manufacturer's brand name (logo), product identification and serial number of the part are given on the back of

the sensor.

	Company logo. MyRay® is a registered trademark of Cefla s.c. - Cefla Dental Group, Imola, Italia
W1234-123	Product identification and serial number of part:

Portable ZEN-X device

The manufacturer's brand name (logo) and product identification are given on the front. Further information is given on the warranty certificate.

	Company logo. MyRay® is a registered trademark of Cefla s.c. - Cefla Dental Group, Imola, Italia
zenX	Name of the appliance

The product identification plate is attached to the back of the interface.


The information contained in the data plate includes:

CEFLA SC - CEFLA DENTAL GROUP VIA BICOCCA 14/C IMOLA (BO) ITALY	Name of manufacturer and production site: Cefla s.c. - cefla Dental Group, Imola, (BO) Italy
	Company logo. MyRay® is a registered trademark of Cefla s.c. - Cefla Dental Group, Imola, Italia
ZEN-X	Name of the appliance
	"Important information" symbol
	"Read the instructions manual" symbol
	"Type B part attached" symbol in accordance with standard EN 60601-1

	ETL mark for compliance to the essential requirements of USA e CANADIAN standards
	The equipment conforms to the requirements stated in EEC directive 93/42 and subsequent changes.
	Symbol "Dispose of carefully as provided for Council Directive 2002/96/CE" (WEEE)
	Class 2 insulation
5V DC USB 500mA	5V supplied directly from Personal Computer USB® 2.0 connection
SN: ZB xxxxx	Part serial number.
 MM/YY	Date of construction (Month / Year+B99)

The following symbols are stamped on the USB® cable:

	USB® 2.0 connection
	Cable in compliance with Standard USB® High Speed® (stamped on cable)

ZEN-Xi device integrated in dental unit

The brand name (logo) of the company that markets the sensor is given on the front. Further information is given on the warranty certificate.

the product nameplate is attached to the edge of the interface.

The information contained in the data plate includes:

	Company logo. MyRay® is a registered trademark of Cefla s.c. - Cefla Dental Group, Imola, Italia
ZEN-Xi	Name of the appliance
ZI xxxxx	Serial number
	The equipment conforms to the requirements stated in EEC directive 93/42 and subsequent changes.

	“Important information” symbol
	“Read the instructions manual” symbol

The operating symbols indicate:

	Stand-by position
	Work position
	Position for taking out the sensor cradle
	Home position

See paragraph “Use of the sensor” for further information.

	The symbol indicating the part can be autoclaved is given on the sensor cradle
--	--

6 Warranty

The manufacturer warrants safety, reliability and performance of the equipment.

 IMPORTANT INFORMATION ! The manufacturer shall not be held liable for any personal injury or property damage arising from failure to heed the following clauses.
--

This warranty is conditional upon the observance of the following prescriptions:

- Heed the conditions given on the warranty certificate found in the box.
- The appliance must be used exclusively in accordance with the instructions contained in this instruction manual.
- The electric installation in the room in which the apparatus is installed must comply to standard IEC 60364-7-710 (Standards regarding electrical installations in medical locations).
- Installation, repairs, settings and in general all the operations for which the equipment covers need to be opened must be carried out only by technicians authorized by the manufacturer.

6.1 Software not covered by warranty

The Software is supplied in its original condition and the Manufacturer shall not be held liable or warrant any original defects or defects originating over time and shall not guarantee quality and proper operation of the software. In addition, the manufacturer shall not honor or provide any warranty regarding conformity of the software to the information given on-line or in electronic documentation or in any case made available except for the warranty for the CD-ROM if damaged or unusable.

Any warranty is also excluded for Software integrated in - or otherwise being a part of - other Software applications developed by third parties. As far as these applications are concerned, the Manufacturer also expressly declares not to have carried out and not to carry out any inspection activity or other activities to guarantee the software operation.

6.2 Limitation of liability.

In no case shall the Manufacturer or its suppliers be responsible for direct or consequential damages (including damages for profit loss or lost earnings or savings, interruption of business activities, loss of data or information or other economic losses) affecting the User or third parties as a consequence of the use or failure to use the Software, also in the event that the Manufacturer had been warned of the possibility of such damages. The present limitation of liability is applicable not just to cases of software use not in compliance with the Manufacturer’s recommendations but also to

cases of software use in compliance with the Manufacturer's recommendations.

7 Minimum system requirements

Before installing the sensor, check for the minimum system requirements listed in the following table:

	PERSONAL COMPUTER: it must be conformed to iec 60601-1 (2nd ed.) regulation or to iec 60950 regulation (information technology equipments safety) and equipped with a insulation transformer.
Component	Requirement
Processor	Minimum: Pentium III ® or similar 800 MHZ Recommended: Pentium IV® 2GHz or higher
Operating System	Windows XP® Home or Professional edition SP 2 - Windows Vista® - Windows 7®
Free hard disk space	At least 150 MB
RAM memory	Minimum: 256 MB Recommended: 512MB or higher
Video board	Video resolution 1024x768 with 16 million colours
Other peripherals	CD-ROM
	USB port ®2.0® HiSpeed®  USB® 1.1 ports should not be used
Screen	Resolution 1024x768 with 16 million colors

X-ray generator	the quality of the images is strictly dependant on the x-ray generator employed. The generator may be AC or DC, power supply 60-70 KV, current 3-10mA. An x-ray tube with focal spot smaller than 0.8mm is recommended. As the sensor is extremely sensitive, a constant potential x-ray generator (DC) with effective time control within tens of milliseconds is recommended. Some older x-ray generators may not be suitable.
Positioner	Use of a positioner guarantees the sensor's sensitive area is centered and that it is perpendicular with the x-ray tube. Use of a positioner is strongly recommended. The clinician should choose the most suitable one based on his/her own experience.
Disposable infection control sheaths	The disposable protection control sheaths are mandatory to safeguard the patient and operator against cross-contamination. The disposable protection control sheaths are class IIa medical devices. They must comply to standard ISO10993-1 regarding biocompatibility and cannot be replaced with others that are poorer quality.

8 Technical characteristics

The technical characteristics and some among the main system functions will largely depend on the characteristics of the radiographic generator and of the software used to display and store the images.

To obtain the best results, it is preferable to use a constant-potential (DC) radiographic generator with long rectangular collimator (focus-to-skin distance not below 30cm).

Old x-ray models that do not permit the exposure times to be sufficiently reduced may not be suitable for use with the ZEN-X / ZEN-Xi system.

Environmental characteristics

- Apparatus designed to work under environmental conditions present in closed work areas in accordance with the parameters set forth by standards I.E.C. 60601-1.
- Not protected against the penetration of liquids (degree of protection (IPX0).
- Never use in the presence of mixtures of flammable anesthetic gas with air, oxygen or nitrogen protoxide.
- Do not install in areas exposed to freezing temperatures.

Operation :

- temperature +10°C - +40°C
- Relative humidity ranging from 20% to 90%
- atmospheric pressure 500 - 1060 hPa

Storage:

- temperature -40 °C - +70 °C
- Relative humidity ranging from 0% to 100%
- atmospheric pressure 500 - 1060 hPa

The apparatus is designed for continuous operation.

General characteristics:

- Class IIa apparatus in accordance with Directive 93/42/CEE and subsequent changes.
- Class 2 apparatus with type B applied part in accordance with standard I.E.C. 60601-1.
- In accordance with European Standards 92/31/CEE and 93/42/CEE and subsequent changes.
- The apparatus complies with European and North American Standards for safety of electrical equipment and electromagnetic interferences therefore dangerous electromagnetic interferences are not generated and the apparatus is not affected by emissions released by equipment during operation.
- Image transfer: through the USB® connection
- Power supply: 5V DC 500mA max, supplied from the PC through the USB® connection
- USB® connection characteristics: USB® 2.0 High Speed ®
- Image transfer: as per TWAIN® standard in multiple image and single mode
- Supported by Windows XP®, Windows Vista®- Windows 7®

ZEN-X interface

- Dimensions: approximately 85 x 54 x 18mm (LxHxT) excluding cables and connectors
- Weight: approximately 50g
- The ZEN-X system interface is not protected against the penetration of liquids (degree of protection IP40).

ZEN-Xi interface

- Dimensions: approximately 162 x 122 x 109 mm (LxHxT) excluding cables and connectors
- Weight: approximately 220g
- The ZEN-Xi system interface is not protected against the penetration of liquids (degree of protection IP20).

9 Sensor characteristics

The sensors feature the following:

- CMOS sensor with Cesium Iodide Scintillator with layer of fiber optics
- Smoothed corners and container with rounded corners for improved ergonomics and patient comfort.
- Cell dimensions 20x20µm
- Theoretical maximum resolution 25 lp/mm
- High sensitivity and dynamic: typical 4mSv saturation, 14 bit converter, 71 db dynamic
- CMOS life: more than 20000 exposures at 60 KV/4 mSv
- CMOS layer of silicon protected with fiber optics against wear caused by direct exposure to X-rays
- Sensor E²PROM memory for storing the operating parameters
- Length of sensor cord: 850 - 1350 mm

Sensor size 1:

- Outside dimensions: 24.9 x 38.9 x 5.3 mm+B102
- Sensitive area: 20 x 30 mm
- Image pixels: 1000 x 1500

Sensor size 2:

- Outside dimensions: 30.4 x 41.9 x 5.7 mm (Width x Height x Thickness at edges)
- Sensitive area: 26 x 34mm
- Image pixels: 1300 x 1700

9.1 Compatibility with x-ray generators

The ZEN-X / ZEN-Xi system works correctly with both conventional x-ray generators called "AC", and newer high frequency generators called "DC". As the sensor is very sensitive, the exposure times need to be reduced compared to

those normally used for conventional x-ray film.

The following chart specifies the exposure times to be used for two different focal distances (according to the collimator used) between the focal spot and skin; the system considers that the distance between the skin and sensor's active surface is approximately 5 cm.

SETTING EXPOSURE TIMES (IN SECONDS)

Length of the cone 12" (30cm)			Length of the cone 8" (20cm)		
					
	0,25	0,16		0,16	0,10
	0,20	0,125		0,125	0,08
	0,16	0,10		0,10	0,063
	0,20	0,125		0,125	0,08
	0,16	0,10		0,10	0,063
	0,20	0,125		0,125	0,08
	0,25	0,16		0,16	0,10

- Tables are referred to a DC generator 60-65 kV, 8 mA, if 4mA is selected it's necessary to double the time.
- If a 70 KV high frequency generator is used, the time given in the table has to be reduced approximately 1/4.
- If edentulous areas are irradiated, the sensor may provide images that are too blackened in the missing areas of the irradiated radiographic subject. In these cases, reduce the time indicated in the table by about 1/4.
- The best results are obtained with a 60 KV high frequency generator, square collimator and 30 cm focus-skin distance (refer to the relevant table).
- To more effectively monitor the distances, the MyRay® centering device equipped with fixed spacer between the centering ring and sensor should be used.
- Before attempting to use the sensor on a patient, practice by taking a few x-ray pictures on unanimated objects with your own x-ray unit.
- Do not exceed the dose in the chart.

MyRay® RX DC HyperSphere, RX DC eXTend

This new generation x-ray generator is specifically designed for use with digital systems.

In fact, it offers 25 different sensitivity settings ("F" factor) to be employed with the digital x-ray receiver used; select the desired setting with the wireless handheld:

- set F=15 for use with the ZEN-X / ZEN-Xi system.

Factor F should be decreased if the images are very dark, according to the positioning technique used.

MyRay® RX AC

This x-ray generator has a selection for digital receivers, shown by a computer icon on the handheld. Make sure the indicator light associated to this function is on. Use the buttons that represent teeth to view the exposure times based on sensitivity factor "S". This factor can be changed with keys + and - on the handheld:

- set mode sensor "on" and S=8 for use with the ZEN-X / ZEN-Xi system.

RX DC - Cefla Dental Group

It is the high frequency x-ray generator prior to MyRay® RX DC HyperSphere. Even though it has only 3 sensitivity factors for traditional film and a setting for digital receivers, it is still possible to save the preferred exposure times.

To correctly use in conjunction with the ZEN-X / ZEN-Xi system, select 60 kV and DIGIT from the control panel. If a 12" long collimator is used, set to 8 mA; with an 8" collimator, set to 4 mA (the light for the selected function should be on).

To save the exposure times suitable for use with the ZEN-X / ZEN-Xi system, select "ADULT patient" and upper molar (18, 17, 16) from the control panel. Use keys + and - to reach the item and set the time to 0.250 s.

Lastly, hold down the save key until the confirmation beep is heard.

10 Disposable covers


IMPORTANT INFORMATION !

The disposable protection control sheaths provided are mandatory to safeguard the patients against cross-contamination. The disposable protection control sheaths are class IIa medical devices. They must comply to standard ISO10993-1 regarding biocompatibility. The disposable protection control sheaths cannot be replaced with others that are poorer quality. Contact the dealer who supplied the sensor to obtain additional disposable protection control sheaths.

Disposable cover application instructions.

- 1) Place a disposable cover complete with its protective sheet on a flat surface. Introduce the sensor through the opening at one end.

- 2) Push the sensor all the way into the disposable cover, paying attention not to break the transparent material.
- 3) If present, remove the protective support sheet.
The operation is now completed.
- 4) After use dispose of the covers as "special" waste.

Use of a positioner guarantees the sensor's sensitive area is centered and that it is perpendicular with the x-ray tube. Use of a positioner is strongly recommended. The clinician should choose the most suitable one based on his/her own experience.

11 Cleaning and disinfecting

Cleaning is the first step of any disinfecting process.

Physically scrubbing with detergents and surface-active substances and rinsing with water removes a considerable amount of micro organisms. If a surface is not clean first, the disinfecting process cannot be successful. If a surface cannot be adequately cleaned, it should be protected with barriers.

IMPORTANT INFORMATION !

The electronic interface's housing is not protected against the penetration of liquids.

Any surfaces that cannot be adequately cleaned should be protected with barriers.

The outside of the apparatus must be cleaned and disinfected with a hospital disinfectant labeled as effective against human immunodeficiency virus (HIV) and hepatitis B virus (HBV) (i.e. low-level) or a hospital tuberculocidal disinfectant (i.e., intermediate-level) specifically designed for small surfaces.

Follow the manufacturer's operating instructions.

The great variety of medical and chemical products used in a dental practice may damage painted surfaces or plastic materials parts. Tests and studies were conducted to prove that surfaces cannot be entirely protected against etching by all the products available commercially. We therefore recommend to use protective barriers whenever this is possible.

The etching effect of chemicals also depends on surface exposure times.

It is therefore important not to leave the product used on the equipment surface to treat for longer than the working time suggested by the manufacturer.

Given the aggressiveness of the active principles used in disinfectants, we recommend to use products containing the following maximum quantities of:

- **Ethanol at 96%.** Concentration: maximum 30 g. for every 100 g. of disinfectant.
- **Propanol.** Concentration: maximum 20 g. for every 100 g. of disinfectant.
- **Combination of ethanol and propanol.** Concentration: combination of the two should be at most 40 g. for every 100 g. of disinfectant.

Cefla sc - Cefla Dental Group performed tests on the compatibility between the most common disinfectants and its own plastic materials.

The test results indicated that the least aggressive agents are:

- ⚠ Incidin Spezial (Henkel Ecolab).
- ⚠ Omnizid (Omnident)
- ⚠ Plastisept (ALPRO) (not tuberculocide as not an alcohol-based disinfectant).
- ⚠ RelyOn Virkosept (DuPont).
- ⚠ Green and Clean SK (Metasys) (not tuberculocide as not an alcohol-based disinfectant)

The same tests showed that the above-mentioned products can therefore be used by adopting, however, the following preATTENTIONS:



IMPORTANT INFORMATION !

- ⚠ Do not use products containing isopropyl alcohol (2-propanol, iso-propanol).
- ⚠ Do not use products that contain sodium hypochlorite (bleach).
- ⚠ Do not use cleaners that contain phenol.
- ⚠ Do not spray the chosen product directly on the equipment surface to disinfect.
- ⚠ Do not mix any of the above products nor mix them with any other products.
- ⚠ Any product should be used in accordance with its manufacturer's instructions.

Cleaning and disinfecting instructions

Clean and disinfect with single-use non-abrasive paper (avoid using recycled paper) or sterile gauze. Sponge cloths or any other re-usable materials are not recommended.

IMPORTANT INFORMATION !

- In order to clean equipment connected to the power mains, shut off the devices and disconnect the power supply

from the outlet before attempting to clean and disinfect the outside.

NOTE: Even if incorporated in the dental chair, the Zen-Xi system works through a USB connection therefore the Personal Computer/WS to which it is connected must be shut off before attempting to remove the sensor cradle.

- Everything used to clean and disinfect must be thrown out after use. Observe current regulations when disposing the material.

Sterilize the sensor cradle only in an autoclave using a sterilization temperature of at least 134°C for 3 minutes.

IMPORTANT INFORMATION !

- **The sensor mount is able to support 500 autoclave sterilization cycles.**
- **Always contact the manufacturer to purchase genuine spare parts.**

The sensor is not designed for autoclave sterilization; use appropriate products to sterilize the outside of it. Use soft disposable paper towel to disinfect the handpiece. Avoid corrosive substances and do not soak in liquids.

IMPORTANT INFORMATION !

The sensor is NOT suitable for autoclave sterilization.

Recommendations for waste disposal

Follow the manufacturer's instructions when disposing whole disinfectant bottles. Do not let the product into the municipal sewer systems and/or waterways.

12 Use of the centering device

The sensor must be held in the correct position using a centering device in order to obtain good x-ray images.

A kit containing centering devices specific for front and rear periapical exposures, bite-wing and endodontics is available in the MyRay® line of products. The individual components in the kit are also available as spare parts. Contact the dealer who supplied the sensor to purchase additional centering devices.

In addition, universal centering devices like RINN® Uni-Grip or KerrHawe® series Bite Sensor or similar can be used. Instructions regarding use of the centering device are given in the Operating instructions included in the centering device kit.

NEVER grasp the sensor with grippers to avoid irreparable damage. Always utilize positioners specifically designed for use with digital x-ray sensors.

13 Installing the sensor

Before installing the ZEN-X / ZEN-Xi system software and drivers, make sure the computer has programs that use TWAIN® for image management (cameras, digital cameras, scanners). Keep in mind that any system drivers installed may interfere with operation of the programs and vice versa. It is recommended to use a dedicated computer for Zen-X. This computer should be used only as a work tool and any software programs that are not needed should be uninstalled.

In order to use the Zen-X or Zen-Xi system the specific sensor software components have to be installed. Consult the iCapture manual and relative instructions.

As regards the ZEN-Xi system integrated in Cefla Dentale dental Group chairs equipped with an incorporated Workstation, not installation procedures are required as all the drivers and software needed are factory installed.

14 Use of the sensor

Power is supplied to the ZEN-X/ ZEN-Xi system directly from the computer's USB® port. Therefore, if a ZEN-X or EN-Xi system is used with an external computer, the computer must be turned on and the electronic interface should be connected to a USB® port. If the ZEN-Xi system is installed in a dental chair equipped with integrated Workstation, the connections have already been made inside the dental chair's system therefore just turn on the workstation.

Insert the type A USB® connector in a free computer USB® port


IMPORTANT INFORMATION !

It is possible to operate only if iCapture is turned on (see the iCapture manual for details about installation and use).


15 Turning the sensor on and off

An x-ray sensor has to be connected to the connector present on the interface in order to use the ZEN-X system.

	Insert the sensor's connector in the interface making sure it is in the right direction. It is correctly inserted if a short click is heard.
---	---

To unplug the connector, squeeze the sides to disengage the lock and remove it avoiding excessive forces or moving sideways.

The x-ray sensor is expensive and fragile, bothered by static electricity. As a result, always remove it from the interface after the interface has been shut off:

- Always use the Stand-By position for ZEN-Xi to disconnect the sensor's connector;
- ZEN-X should be shut off with the specific key or by disconnecting it from the USB® port.

16 Acquiring an x-ray image

An x-ray unit has different capturing modes depending on the acquisition software used. Consult the instructions supplied by the manufacturer of the software used.

If the iCapture Monitor software supplied with the ZEN-X system is employed, use is extremely simple.

Boot an image acquisition program by selecting image acquisition from iCapture.

	Do not take x-ray pictures on a patient when testing the system for the first time or verifying correct operation. Use phantoms to conduct tests.
--	--


Put the x-ray sensor in the patient's mouth.

Make sure the status light is green.

Adjust the exposure times

Take X-ray picture


The image will appear on the computer screen after a few seconds, and, if enabled, it can be seen in the iCapture Monitor preview window and right-hand column next to the main window.

No other steps are required to capture further images after the first x-ray picture has been taken.

The last image captured is shown in the preview window.

The new images will appear under the first one in the right-hand column of the main window.

17 Status light

Always check system status before attempting to take x-rays on a patient. Make sure the operating light is GREEN before taking an x-ray on a patient.

A red light indicates a serious fault for which it is not possible to take the x-ray (the self-diagnosis system detected at start up).

A yellow light indicates a temporary situation (initialization, image transfer) or a fault which the user should be aware of (e.g.: sensor not connected);

A blue light indicates the system is connected to the Personal Computer but not on;

Consult the sensor status table for a detailed description of the messages.

SENSOR STATUS TABLE

Red	Error during execution of self-diagnostic testing
Yellow flashing light	• Sensor not connected • iCapture not on
Steady yellow light	• Setup in progress • Calibration in progress
Blue (slow flashing light)	Sensor in stand-by
Green	Ready: sensor on and ready to receive an image

18 ZEN-X

ZEN-X does not have its own start button. Once connected to a PC, the TWAIN® driver is automatically started which will initialize the ZEN-X system if iCapture Monitor is running (factory settings).

Initially a yellow light will come on and “WAIT” will appear on a yellow background in the MyRay® Driver TWAIN® window in ZEN-X.

If the sensor is connected, “WAIT” appears in yellow in window MyRay® Driver TWAIN® and after a few seconds “READY” appears in green in window MyRay® Driver TWAIN® and the light on ZEN-X will become green.

At this point, the ZEN-X system is ready to receive an x-ray image.

If the sensor is NOT connected, “SENSOR NOT CONNECTED” appears in red in window MyRay® Driver TWAIN® after a few seconds and the light on ZEN-X will become blue (stand-by).

If the ZEN-X system is connected to a PC in which MyRay® Driver TWAIN® is NOT on (the factory settings have been modified by the user), it is not possible to acknowledge and initialize. In this case, the light will become blue to signal the ZEN-X system is in stand-by and x-ray pictures cannot be taken.

Consult the iCapture manual to start MyRay® Driver TWAIN®.

When MyRay® Driver TWAIN® is turned on, the drivers are automatically identified and started and the system is ready to function after a few seconds.

19 ZEN-Xi

As ZEN-Xi is installed on the instrument board as the sixth instrument, it may be to the right or left of the instrument board depending in the particular case. For this reason, two indicator lights are provided. They are placed so that the operator can always see at least one of them under regular work conditions.

The sensor cradle can be turned and placed in four distinct positions.

Operation is similar to that of the ZEN-X system when in this position.

Once connected to a computer with MyRay® Driver TWAIN® started, the required drivers are automatically identified and initialized.

Initially a yellow light will come on and “WAIT” will appear on a yellow background in the MyRay® Driver TWAIN® win-

dow in ZEN-Xi.

If the sensor is connected, "READY" appears in green in window MyRay® Driver TWAIN® after a few seconds and the light on ZEN-X will become green.

At this point, the ZEN-X system is ready to receive an x-ray image.

If the sensor is NOT connected, "SENSOR NOT CONNECTED" appears in red in window iCapture Monitor after a few seconds and the light on ZEN-X will become yellow.

To quit this condition, connect a sensor to the interface or switch over to position Door closed.

Stand-By

This position switches ZEN-Xi to Stand-by.

The light on ZEN-Xi becomes blue and the system switches to stand-by. When turned to On, the system is immediately read to operate.

	Always use the Stand-By position to disconnect the sensor's connector.
---	--

Switching from Stand-by to On is detected by the electronic interface that automatically opens the connection between the sensor and computer or workstation.

Door closed

As ZEN-Xi is physically attached to the dental chair instrument board, dirt and sprays may hit the connector when the sensor is not connected to the interface.


	Turn the sensor cradle until the opening is covered by the relative part for hygienic reasons and to avoid damaging electronic parts.
---	---

If the door is left opened without the sensor, the light will become yellow and flash and "sensor not connected" will appear in window iCapture Monitor in red.

Extraction

The sensor cradle can be taken out of the interface in this position so that cleaning and disinfecting operations can be performed.

Clean and disinfect the sensor cradle on a regular basis. The sensor mount can be sterilized in autoclave. It is recommended to have an adequate number of sensor cradles so that a clean one can be used for each patient for whom x-rays need to be taken. Additional cradles can be ordered from the dealer who supplied the apparatus.

20 Quality of the x-ray images

Unlike common x-ray film, the ZEN-X / ZEN-Xi system tends to automatically correct any exposure errors, providing images that are always usable. Notwithstanding the sensor allows images with a wide range of grey levels to be captured, computer monitors display only 256, therefore in most cases the software will obtain a satisfactory image even from a picture taken with incorrect exposure. Keep in mind there are however boundaries beyond which results cannot be corrected.

The x-ray sensor is more sensitive than x-ray film therefore exposure times usually have to be reduced. See the information given in paragraph **Sensor characteristics - Compatibility with x-ray generators**.

To obtain top performance of digital x-ray sensors, it is important to keep in mind that there are some differences compared to film. Insufficient exposure is clearly seen on x-ray film as the areas corresponding to soft tissues are dark. On the other hand, when a digital sensor is used, the image's background noise increases (salt and pepper noise) and the tonal range is insufficient.

Excessive exposure (time too long) on x-ray film causes the image to be too dense (dark) while the image will lose contrast in a digital sensor.

It is a common mistake to confuse excessive exposure with insufficient exposure therefore further increase the exposure times.

It is important to check and note this limit with your own x-ray system in order to be sure not to exceed it during dental treatment as the quality of the images obtained under these conditions will be poor or even unusable.

IMPORTANT INFORMATION !

Before attempting to take x-ray pictures on patients, it is advisable to take some test pictures on phantoms, comparing the results obtained to the usual ones. Identify the best exposure conditions for your own x-ray system through trial and error.

21 Problem solving

Problem found	Possible causes	Remedies
Doubts concerning the sensor efficiency	Falling, banging, malfunctioning	Do not use the sensor on a patient. Conduct tests by capturing x-ray images using a phantom. If there are doubts about proper operation, do not use the sensor and contact the Service Center.
Image loss	Error in the management program or PC operating system	The last image captured can be recovered by downloading it again from the sensor from window "TWAIN® data source" (see paragraph Recovering the last image captured iCapture user instructions). Do not shut off the PC or disconnect the interface from the USB® port to avoid permanently losing the image.

The system does not come on.	USB® cable not connected	Connect the USB® cable to a port.
The system does not come on.	Personal Computer USB® cable or USB® port defective	Check the personal Computer's USB® cable and USB® port with another device, e.g. a pen drive. test the equipment on another PC. The status light should be on (flashing yellow) even without installing the software.
The system does not come on.	ZEN-X fault	Do not use the sensor, contact the Technical Service Center.
The status light on the interface stays red	ZEN-X fault or sensor damaged or defective	Do not use the sensor, contact the Technical Service Center.
Sensor does not start; the light is yellow and flashes	Driver not present, defective or damaged. New software has been recently installed in the PC or the PC has been used with external connections (Internet).	Run an antivirus check. Reinstall the software. Only use your PC as a working tool, avoiding connection to external networks.
ZEN-Xi cannot be started (it stays in Stand-By) the light is flashing yellow	The sensor cradle is not correctly inserted or is not turned to the correct position.	Check rotation and verify the sensor cradle is in the correct position.
The system comes on but the light stays yellow and flashes: a fault message appears on the computer.	Quality of USB cable not good enough or cable is too long. NOTE: the maximum length of a top quality USB cable is about 4.5 m.	Replace the USB cable; eliminate any extensions; try to use a HUB energized in the last section before the x-ray system connection.
The Zen-X x-ray system is not acknowledged	A version of iCapture released prior to the introduction of the Zen-X x-ray system is installed in the PC.	The Zen-X system is found starting from iCapture version 2.2. Close iCapture and install an updated version.
An ERROR message followed by a number appears on the PC	Sensor or interface malfunction	Write down the message and inform technical personnel. Do not use the sensor, contact the Technical Service Center.
An image is acquired but with poor tone range and/or massive background noise.	Underexposed image.	Use a longer exposure time, make sure the x-ray generator works correctly.
The image is captured but it is very light with little contrast.	Over-exposed image	Use a shorter exposure time, check the settings of the x-ray generator.

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